

# Data Transmission

To ensure reliable, uninterrupted transmission of data, you should pay attention to the following:

## Receiver Out of Range

- If the receiver is too far from the Bravo capsule, a beep will be heard for 30 seconds and the display will flash C1 or C2 (P1 or P2 may flash depending on receiver version). This indicates that the receiver is out of range and may not be receiving information from the pH capsule. Move the receiver to your chest until the beep stops and the C1 or C2 is no longer visible.
- If the C1 or C2 display continues or if you have any problems or questions during your test, call your doctor.
- At the completion of your pH test, return the Bravo pH receiver and your diary to **Huron Gastro** on...

**NOTE:** Stay within one meter (3 feet) of the receiver during your pH test!

## External Interference

As with any radiotelemetry device, the Bravo pH capsules signal may be interrupted by other electronic or electrical devices operating at the same or similar frequencies.

**Solution:** If the Bravo receiver beeps and the message C1 or C2 appears on the display, you should move away from device(s) that appear(s) to cause the beep and C1/C2 message.

# Patient Diary

During the test, it is important that you eat, drink, work, and exercise as you normally would so your doctor gets reliable data to diagnose your condition. Record meals, sleep, and other important activities during your test in your Patient Diary using the time on the receiver's display.




## Using the Patient Diary

During your pH test, you should record your symptoms (Heartburn, Regurgitation, Chest Pain, Coughing, etc.) in the Patient Diary.

## Using the Symptom Buttons

Besides using the Patient Diary, you may also record your symptoms electronically on the receiver by pressing the appropriate symptom button.

A green indicator light remains lit for 3 seconds after a button is pushed. You may hear a beep when a symptom button is pushed. Not all clinicians select this feature: therefore your receiver may not beep.

Symptom Button	Symptom
	Heartburn
	Regurgitation
	Chest Pain